

1

2

3

4

A

**PLASTER TROLLEY - SERIES 300**

ISO MODULES 600 x 400 mm  
 POLYETHYLENE WORK SURFACE  
 WITHSTANDS ALCOHOL, CHLORHEXIDINE  
 (SEE WEB - TECHNICAL FEATURES)

**49RI** - 2 BASINS 9 LITRES AT RAIL

- AISI-304 STAINLESS STEEL SUPPORT,  
REMOVABLE FOR CLEANING
- 2 POLYPROPYLENE BINS, YELLOW COLOUR

**CARRO DE YESOS - SERIE 300**

MÓDULOS ISO 600 x 400 mm  
 SUPERFICIE DE TRABAJO DE POLIETILENO  
 RESISTE ALCOHOL, CLORHEXIDINA  
 (VER WEB - CARACTERÍSTICAS TÉCNICAS)



**49RI** - 2 CUBETAS 9 LITROS A RIEL

- SOPORTE ACERO INOXIDABLE AISI-304,  
DESMONTABLE PARA LIMPIEZA
- 2 CUBETAS POLIPROPILENO, COLOR AMARILLO

B

C

D




640 x 480 mm

333 W.26(R2.R2).30.37L.49RI

900 x 630 mm

N322 W.26(R2N).37NR.49RI

06		Dibujado	Fecha	ESCALA	
05		OT1	09-06-17		
04		Comprobado	Fecha		
03					
02					PLANO Nº
01					991-0282-00
MODIFICACIONES		MOD. 333 W.26(R2.R2).30.37L.49RI			