

TECHNICAL FEATURES

SIMPLE TROLLEY 100 SERIES

640 x 480 mm

EXTENSIVE MOBILITY TROLLEYS -HOSPITALIZATION UNITS

WORK SURFACE

- Polyethylene, PE, high resistance: 100% pure alcohol (not diluted), bleach, detergent, soap solutions, oils, polish, aqueous ammonia, chlorhexidine.
- Rounded shapes, asepsis, hygiene, easy-to-clean.
- Protection on all 4 sides, including the front.

<http://www.insaustimedicaltrolleys.com/en/technical-features/work-surface/>



MONOBLOC STRUCTURE

- AISI-304 stainless steel.
- Fine satin finish, asepsis, hygiene.
- Its compact, sturdy structure.

<http://www.insaustimedicaltrolleys.com/en/technical-features/monobloc-structure/>

ISO DRAWERS 600 x 400 mm

- ABS ISO drawers of standard dimensions according to ISO 3394, compatible and interchangeable with ISO furniture of hospitals using this system
- Heights 50, 100, 200 mm.
- The drawers supported by a ISO drawer framework.
- Total opening, steel ball-bearing slides with self-closing system
- Front stop, it protects the slides against strong impacts when closing.
- Card holder with identification cards ISO norm 7810.

<http://www.insaustimedicaltrolleys.com/en/technical-features/drawers/>



Ergonomic
push handle

Polyamide handle,
PA, ergonomic
Easy-to-clean

Polyethylene drawer
front part, PE
6 colours.

WHEELS

Ø 125 mm

2 with brake
Ergonomic pedal

2 Antistatic

1005

