

## **EMERGENCY - SIMPLE TROLLEY 100 SERIES** **1299 R**

**DIMENSIONS** 640 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 1005 mm (work surface – floor).

### **WORK SURFACE**

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

### **MONOBLOC STRUCTURE**

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

### **ISO DRAWER FRAMEWORK**

- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Interior frame ISO 600 x 400 mm.
- Polyethylene drawer front part, thickness 10 mm, 6 colours (Ref. G, Y, W, B, F, R).
- Ergonomic handle.
- Card holder with white cards, removable.

### **ISO DRAWERS 600 x 400 mm**

- Heights 50, 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.

### **EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION**

- Elevated push handle. Ergonomics, height 1045 mm.
- Double wheels Ø 125 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at push handle side.



## **1299 R**

**165 R.12.14.28.32R.34.36.39R.40**

### **ACCESSORIES**

**12** MANUAL LOCKING

**14** NUMBERED DISPOSABLES SEALS - 100 units

**28** ROTATING ARM FOR DEFIBRILLATOR

**32R** SIDE SUPPORT

**34** CARDIAC MASSAGE BOARD

**36** GAS BOTTLE STAND

**39R** SUCTION PUMP SUPPORT

**40** SERUM HOLDER