

## **ANAESTHESIA - SIMPLE TROLLEY 300 SERIES**

### **3685 B**

**DIMENSIONS** 640 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 980 mm (work surface – floor).

#### **WORK SURFACE**

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.
- Raised sides make it possible to move the trolley without needing the push handle.

#### **MONOBLOC STRUCTURE**

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

#### **ISO DRAWER FRAMEWORK**

- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Interior frame ISO 600 x 400 mm.
- Polyethylene drawer front part, thickness 10 mm, 6 colours (Ref. G, Y, W, B, F, R).
- Ergonomic handle.
- Card holder with white cards, removable.

#### **ISO DRAWERS 600 x 400 mm**

- Heights 50, 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.

#### **STATIC TROLLEYS**

- Without push handle - Easy left panel cleaning.
  - Free left side for mounting accessories.
- Double wheels Ø 100 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at the front part.



### **3685 B**

**365 B.26(RA4.A5.R2).30.(2)32L**

#### **ACCESSORIES**

**26(RA4.A5.R2)** BINS DISPENSERS - CENTRAL SHELF

**30** FOLD DOWN TRAY

**(2)32L** SIDE SUPPORTS