

## **BROSELOW PAEDIATRIC EMERGENCY - SIMPLE TROLLEY 100 SERIES** **1200 W**

**DIMENSIONS** 640 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 1005 mm (work surface – floor).

### **WORK SURFACE**

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

### **MONOBLOC STRUCTURE**

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

### **ISO DRAWER FRAMEWORK**

- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Interior frame ISO 600 x 400 mm.
- Polyethylene drawer front part, thickness 10 mm.
- Ergonomic handle.

### **ISO DRAWERS 600 x 400 mm**

- 6 drawers 100 mm heights. Interchangeable, ABS.
- ISO norm 3394.
- Upper drawer without divider.
- 4 drawers divided into 2 equal parts (1 Ref. DS10 per drawer).
- Auxiliary lower drawer without divider.

### **EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION**

- Elevated push handle. Ergonomics, height 1045 mm.
- Double wheels Ø 125 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at push handle side.



### **1200 W**

**164M6 W.HBR.12B.14.36.40.(4)DS10**

### **ACCESSORIES**

**HBR** BROSELOW 9 COLOUR CARD HOLDER SET

**12B** INDIVIDUAL MANUAL LOCKING

**14** NUMBERED DISPOSABLE SEALS - 100 units

**36** GAS BOTTLE STAND

**40** SERUM HOLDER

**(4)DS10** ISO DRAWER DIVIDER SETS