

1

2

3

4

**SUTURE TROLLEY****SERIES 300**

DIMENSIONS 640 x 480 mm

ISO MODULES 600 x 400 mm

POLYETHYLENE WORK SURFACE

WITHSTANDS ALCOHOL, CHLORHEXIDINE  
(SEE WEB - TECHNICAL FEATURES)**R2** - CENTRAL SHELF, VARIABLE HEIGHT  
OPTIONAL: 3 OR 4 SHELVES FOR  
SUTURES BOXES**CARRO DE SUTURAS****SERIE 300**

DIMENSIONES 640x 480 mm

MÓDULOS ISO 600 x 400 mm



SUPERFICIE DE TRABAJO DE POLIETILENO  
RESISTE ALCOHOL, CLORHEXIDINA**R2** - ESTANTE CENTRAL, ALTURA VARIABLE  
OPCIONAL: 3 O 4 ESTANTES PARA  
CAJAS DE SUTURAS

A

B

C

D



06		Dibujado	Fecha	ESCALA	
05		OT1	02-08-19		
04		Comprobado	Fecha		
03					
02					PLANO Nº
01					991-0242-00
MODIFICACIONES		MOD. 365 W.26HE(R2.R2.R2.R2)			