

1

2

3

4

SUTURE TROLLEY**SERIES 300**

DIMENSIONS 900 x 630 mm

ISO MODULES 600 x 400 mm

POLYETHYLENE WORK SURFACE

WITHSTANDS ALCOHOL, CHLORHEXIDINE
(SEE WEB - TECHNICAL FEATURES)**R2N** - CENTRAL SHELF, VARIABLE HEIGHT
OPTIONAL: 3 OR 4 SHELVES FOR
SUTURES BOXES**CARRO DE SUTURAS****SERIE 300**

DIMENSIONES 900 x 630 mm

MÓDULOS ISO 600 x 400 mm




SUPERFICIE DE TRABAJO DE POLIETILENO
RESISTE ALCOHOL, CLORHEXIDINA**R2N** - ESTANTE CENTRAL, ALTURA VARIABLE
OPCIONAL: 3 O 4 ESTANTES PARA
CAJAS DE SUTURAS

A

B

C

D



06		Dibujado	Fecha	ESCALA	
05		OT1	01-08-19		
04		Comprobado	Fecha		
03					
02		MOD. N365 W.26HE(R2N.R2N.R2N.R2N)			PLANO Nº
01					991-0243-00
MODIFICACIONES					