

TREATMENT - SQM SQUARE TROLLEY Q158 W

DIMENSIONS 480 x 480 mm

ADEQUATE HEIGHT, ERGONOMICS 1005 mm (work surface – floor).

WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

SQUARE DRAWERS

- Polypropylene, PP. 3 colours - B Blue, W White, R Red.
- Heights 57, 130 mm.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 10 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION

- Elevated push handle. Ergonomics, height 1045 mm.
- Double wheels Ø 125 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at push handle side.



Q158 W

SQ161 W.26(RA5Q.A5.R1Q).30.33L.46.50

ACCESSORIES

26(RA5Q.A5.R1Q) BINS DISPENSERS - CENTRAL SUPPORT

30 FOLD DOWN TRAY

33L DISPOSABLE CONTAINER SUPPORT

46 3 GLOVES BOX HOLDER

50 13 LITRES BIN