

## **EMERGENCY - SIMPLE TROLLEY 100 SERIES**

### **1224 R**

**DIMENSIONS** 640 x 480 mm  
**ADEQUATE HEIGHT, ERGONOMICS** 1005 mm (work surface – floor).

#### **WORK SURFACE**

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

#### **MONOBLOC STRUCTURE**

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

#### **ISO DRAWER FRAMEWORK**

- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Interior frame ISO 600 x 400 mm.
- Polyethylene drawer front part, thickness 10 mm, 6 colours (Ref. G, Y, W, B, F, R).
- Ergonomic handle.
- Card holder with white cards, removable.

#### **ISO DRAWERS 600 x 400 mm**

- Heights 50, 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.

#### **EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION**

- Elevated push handle. Ergonomics, height 1045 mm.
- Double wheels Ø 125 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at push handle side.



### **1224 R**

165 R.12.14.26LW(R7).28EV.32R.34.36.40

#### **ACCESSORIES**

**12** MANUAL LOCKING

**14** NUMBERED DISPOSABLES SEALS - 100 units

**26LW(R7)** CENTRAL RAIL

**28EV** ROTATING ARM

**32R** SIDE SUPPORT

**34** CARDIAC MASSAGE BOARD

**36** GAS BOTTLE STAND

**40** SERUM HOLDER