

## MULTIFUNCTION - SQS SQUARE TROLLEY Q077 W

**DIMENSIONS** 480 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 980 mm (work surface – floor).

### WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.
- Raised sides make it possible to move the trolley without needing the push handle.

### MONOBLOC STRUCTURE

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

### SQUARE DRAWERS

- Polypropylene, PP. 3 colours - B Blue, W White, R Red.
- Heights 57, 130 mm.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 10 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

### STATIC TROLLEYS

- Without push handle - Easy left panel cleaning.
  - Free left side for mounting accessories.
- Double wheels Ø 100 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at the front part.



### Q077 W

SQ061 W.26LW(R1Q).30.33L

### ACCESSORIES

**26LW(R1Q)** CENTRAL SUPPORT

**30** FOLD DOWN TRAY

**33L** DISPOSABLE CONTAINER SUPPORT