

TRANSPORT TROLLEYS - TROLLEY WITH ISO PANELS 620 DN.(12)72.(4)73

DIMENSIONS 900 x 630 mm

HEIGHT 1645 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Union with rivets and TIG welding.

ISO PANELS

- Flexible organisation inside the trolley.

ISO DRAWERS 400 x 600 mm

- Heights 50, 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.
- Maximum capacity 10 + 10 = 20 ISO drawers 100 mm height.




2 DOORS

- Total opening 270°.
- Phenolic, 10 mm thickness.
- Lock with 4 keys.
- Card holder with white cards, removable, left door.
- DN Light grey.

VERTICAL HANDLES

- 2 left side, 2 right side.

EXTENSIVE MOBILITY TROLLEYS

- 2 bumpers for hits and 4 wheel bumpers.
 - Ø 125 mm wheels.
 - 2 antistatic, 2 with dual lock at the right side.
- Pedal position:  Fixed wheel.  Swivel wheel.  Locked wheel.



620 DN.(12)72.(4)73

ACCESSORIES

(12)72 ISO DRAWERS 100 mm

(4)73 ISO DRAWERS 200 mm