

## TRANSPORT TROLLEYS - TROLLEY WITH SLIDES 657 DN

**DIMENSIONS** 640 x 480 mm  
**HEIGHT** 1645 mm

### MONOBLOC STRUCTURE

- Stainless steel AISI-304. Union with rivets and TIG welding.

### ISO DRAWER FRAMEWORK

- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Interior frame ISO 600 x 400 mm.

### ISO DRAWERS 600 x 400 mm

- Heights 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.




### 2 DOORS

- Total opening 270°.
- Phenolic, 10 mm thickness.
- Lock with 4 keys.
- Card holder with white cards, removable, left door.
- DN Light grey.

### VERTICAL HANDLES

- 2 left side, 2 right side.

### EXTENSIVE MOBILITY TROLLEYS

- Ø 125 mm wheels.
  - 2 antistatic, 2 with dual lock at the right side.
- Pedal position:  Fixed wheel.  Swivel wheel.  Locked wheel.

