

TECHNICAL FEATURES

DOUBLE TROLLEY 300 SERIES

900 x 630 mm

STATIC TROLLEYS - EMERGENCY ROOM, OPERATING ROOM, REANIMATION AREA, ICU

WORK SURFACE

- Polyethylene, PE, high resistance: 100% pure alcohol (not diluted), bleach, detergent, soap solutions, oils, polish, aqueous ammonia, chlorhexidine.
- Protection on all 4 sides, including the front.

<http://www.insaustimedicaltrolleys.com/en/technical-features/work-surface/>



ISO DRAWERS 400 x 600 mm

- ABS ISO drawers of standard dimensions according to ISO 3394, compatible and interchangeable with ISO furniture of hospitals using this system
- Heights 50, 100, 200 mm.
- The drawers supported by a ISO drawer framework.
- Total opening, steel ball-bearing slides with retention system
- Front stop, it protects the slides against strong impacts when closing.

<http://www.insaustimedicaltrolleys.com/en/technical-features/drawers/>



Trolley prepared to mount on accessories at destination

Monobloc structure of AISI-304 stainless steel

Polyamide handle, PA, ergonomic
Easy-to-clean

Polyethylene drawer front part, PE
6 colours

WHEELS
Ø 100 mm

2 Antistatic

2 with brake
Ergonomic pedal

