

TECHNICAL FEATURES

SQS SQUARE TROLLEY

480 x 480 mm

STATIC TROLLEYS - EMERGENCY ROOM, OPERATING ROOM, REANIMATION AREA, ICU

WORK SURFACE

- Polyethylene, PE, high resistance: 100% pure alcohol (not diluted), bleach, detergent, soap solutions, oils, polish, aqueous ammonia, chlorhexidine.
- Protection on all 4 sides, including the front.

<http://www.insaustimedicaltrolleys.com/en/technical-features/work-surface/>



SQUARE DRAWER

- Polypropylene drawers, PP.
- Resistance to chemicals similar to the polyethylene.
- Colour white, blue, red.
- Total opening, steel ball-bearing slides with retention system.
- Front stop, it protects the slides against strong impacts when closing
- Identification cards ISO norm 7810.

<http://www.insaustimedicaltrolleys.com/en/technical-features/drawers/>



Trolley prepared to mount on accessories at destination

Monobloc structure of AISI-304 stainless steel

Polyamide handle, PA, ergonomic
Easy-to-clean

WHEELS
Ø 100 mm

2 Antistatic

2 with brake
Ergonomic pedal

