

1

2

3

4

STORAGE FURNITURE**700 SERIES**

DIMENSIONS: 900 x 630 mm

HEIGHT: 1600 mm

MONOBLOC STRUCTURE OF
AISI-304 STAINLESS STEEL

PHENOLIC DOORS WITH KEY LOCKING

62 - ADJUSTABLE SHELF 400 x 610 mm**72** - ISO DRAWER 100 mm 400 x 600 mm**73** - ISO DRAWER 200 mm 400 x 600 mm

MAX.CAPACITY 10 + 10 = 20 ISO DRAWERS 100 mm

MOBILIARIO DE ALMACÉN**SERIE 700**

DIMENSIONES 900 x 630 mm


ALTURA: 1600 mm

ESTRUCTURA MONOBLOQUE DE
ACERO INOXIDABLE AISI-304


PUERTAS DE FENÓLICO CON BLOQUEO DE LLAVE

62 - BANDEJA REGULABLE 400 x 610 mm**72** - CAJÓN ISO 400 x 600 mm DE 100 mm**73** - CAJÓN ISO 400 x 600 mm DE 200 mm

CAPACIDAD MÁX. 10 + 10 = 20 CAJONES ISO DE 100 mm



975

06		Dibujado	Fecha	ESCALA	
05		OT1	23-03-22		
04		Comprobado	Fecha		
03					
02		MOD. 720 DN.(3)62.(6)72.(2)73		PLANO Nº	
01				991-0725-00	
MODIFICACIONES					